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Regulations and Rules for use of Fluorescence Microscopes in LSII-332

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if you are already a qualified user of one fluorescence microscope

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and you want to use another one, you need to:

6. 由合格使用者全程陪同，自己獨立進行三次實驗，總時數達3小時以上(陪同之合格使用者需簽名)。
operate with a certified user aside three times, for a total of 3 hours (signatures of a certified user are required)
7. 請指導老師寄發更新的貴重儀器使用同意書。
email TA an updated signed Agreement Sheet of Instrument by the lab PI.
8. 尋找管理者核對使用記錄與進行操作測驗。使用紀錄、操作測驗與同意書均核對無誤後，即成為另一台顯微鏡的合格使用者。
check all the records with TA and take an operation exam. After everything is done, the user will be qualified for another microscope.

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Before using the microscopes, new users are encouraged to attend the microscopy training workshop or have taken any related course about microscopy (e.g., Prof. DT Chang's biological macromolecules, etc.) to ensure a full understanding of the basics of microscopy,

II. 操作不當之處罰 Punishment about misconduct

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Users must book time block each time before use. Users must also fill the form in the notebook after using the microscope. Users have to cancel the reservation if users cannot use the microscope on time.

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Users without certification are not allowed to use reservation system.

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Users will be subjected to one offense for disobeying the above rules.

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The users should also fill the form in the notebook even if using the computer only.

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Backing up data in a separate storage device immediately is strongly encouraged.

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Files will be removed from the hard drive periodically at the beginning of each month, and will be announced in advance.

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If you find any abnormality in operating the microscope, you need to

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In case of emergency, call the administrator (LY Lin) immediately.

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Users are responsible for the cleaning of the microscope room, including the floor, the desk, the external part of the microscope, and the oil objectives, etc..

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Do not eat or drink in the microscope room.

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1. 各顯微鏡之維修與消耗性費用(請見附錄)將由使用實驗室按時間比例分攤。
The maintenance fee and cost of consumable parts of each microscope (please see appendix) will be partitioned to each laboratory according to total using time.
2. 帳單會以電子郵件寄到各指導教授及院秘書，費用將直接由教授的系所經費中扣除。
The bill will be emailed to both individual PI and secretary, and money will be deducted from PI's departmental account directly.

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Basic cost of consumable items and maintenance

Consumable part	Unit price	
Filter set (about 3 year; about 8000 hr)	24000	Upright: 4 sets Inverted: 3 sets
Filter module (about 3 year; about 8000 hr)	80000	Stereo: 2 sets;
Hg lamp 100W (250 hr)	10000	One per microscope
Regular maintenance (once per year)	10000	Yearly per microscope